

# MagCore® Nucleic Acid Extraction Kits



# MagCore® Automated Nucleic Acid Kits Selection Guide



	Cat.No.	Cat.No.	Cat.No.	Cat.No.	Whole Blood	Buffy Coat	Plasma/Serum	Cell-free body fluids	Urine	Free circulating DNA	Amniotic Fluid	Cultured Cells	Animal Tissue	Plant Tissue	Bacteria/Sputum	Swab	FFPE	Stool	Forensic-Specimens
Genomic DNA	<b>101</b> MagCore® Genomic DNA Whole Blood Kit (Speedy installation)	MGB400-01	MGB400-02		■	■													
	<b>102</b> MagCore® Genomic DNA Whole Blood Kit	MGB400-03	MGB400-04		■	■													
	<b>104</b> MagCore® Genomic DNA Large Volume Whole Blood Kit (1.2 ml)		MGB1200		■														
	<b>105</b> MagCore® Plasma DNA Extraction Kit (1.2 ml)		MPD1200				■	■	■										
	<b>115</b> MagCore® Circulating DNA large volume kit (3-4 ml)	MPD4000-01		MPD4000-03			■	■	■										
	<b>106</b> MagCore® Genomic DNA Whole Blood Kit (For Genotyping)	MGB400-07	MGB400-08		■														
	<b>110</b> MagCore® Cultured Cells DNA Kit	MCC-01	MCC-02							■	■								
	<b>301</b> MagCore® Genomic DNA Plant Kit	MGP-01	MGP-02											■					
	<b>401</b> MagCore® Genomic DNA Tissue Kit	MGT-01	MGT-02											■		■	■	■	■
	<b>405</b> MagCore® Genomic DNA FFPE One-Step Kit	MGF-01	MGF-03														■		
	<b>406</b> MagCore® Forensic DNA Direct Kit		MFC-03																■
<b>502</b> MagCore® Genomic DNA Bacterial Kit	MBB-01	MBB-02													■				
Viral Nucleic Acids	<b>201</b> MagCore® Viral Nucleic Acid Extraction Kit	MVN400-01	MVN400-02				■	■	■									■	
	<b>202</b> MagCore® Viral Nucleic Acid Extraction Kit (Low PCR Inhibition)	MVN400-03	MVN400-04				■	■	■									■	
	<b>203</b> MagCore® Viral nucleic acid extraction kit (high sensitivity)	MVN400-05	MVN400-06				■	■	■									■	
	<b>210</b> MagCore® Viral Nucleic Acid Large Volume Extraction Kit (2.4 ml)		MVN2400				■	■	■									■	
	<b>211</b> MagCore® Viral Nucleic Acid Large Volume Extraction Kit (1.2 ml)		MVN1200				■	■	■									■	
Total RNA	<b>601</b> MagCore® Total RNA Whole Blood Kit	MRN-01	MRN-02		■														
	<b>605</b> MagCore® Total RNA FFPE One-step Kit	MRF-01	MRF-03															■	
	<b>610</b> MagCore® Total RNA Cultured Cells Kit		MRC-01	MRC-02								■							
	<b>631</b> MagCore® triXact RNA Kit	MRX-01	MRX-03			■						■	■	■					